

SuperLine & Implantium

Surgical / Prosthetic Manual

DentiumUSA
Developed by Clinicians for Clinicians

SuperLine[™] & IMPLANTIUM[®]

*A New Choice
For the Customer*





SURGICAL MANUAL

Surgical Drill Sequence	04
Drilling Depth Guide	06
Fixture Connection	08
Installation Procedure	09
Surgical Kit Maintenance	09

Surgical Drill Sequence


Drilling Sequence Guide



During fixture insertion, 30~45N·cm torque at 20rpm is recommended.

- Countersink drill is used in cases with dense cortical bone.
- The countersink drill's actual diameter is 0.1mm larger than the fixture platform.

Determination of Fixture Top Level




- It is recommended that the top level of the fixture be located 0.5mm below the crestal bone.

Depth Indication



- Use the Depth Gauge after first drill / Lindemann first drill to check depth of drilling
- Place the Depth gauge against the wall of the osteotomy

Drilling Depth Guide




Countersink Depth Guide

- Drilling Depth of the countersink depends on the patient's bone quality.
- If the bone density is D1~D2, it is recommended to drill up to the top line (I) of laser mark on the countersink.
- If the bone density is D3~D4, it is recommended to drill up to the bottom line (II) of laser mark on the countersink.


Platform: $\varnothing 3.6$ / Body: $\varnothing 3.4$

(1000rpm/30~45N·cm)




Platform: $\varnothing 4.0$ / Body: $\varnothing 3.8$




(1000rpm/30~45N·cm)



Platform: $\varnothing 4.5$ / Body: $\varnothing 4.3$

(1000rpm/30~45N·cm)





※ Note : Instead of Lindemann guide drill, Guide drill may be used
 Instead of Lindemann first drill, First & Pilot drill may be used.
 * Do not use Pilot drill when placing 3.6 platform / 3.4 body fixture.


Fixture Connection



Caution_ When opening the fixture pack, hold the fixture container upward and engage the adapter into the fixture.




By handpiece
20rpm/35N · cm



By ratchet



Directions Using the Hand-piece / Ratchet Adapter



Hand-piece adapter

Ratchet adapter

The adapters and the fixture internal hex must be connected firmly together.

Installation Procedure

Cover Screw




By hex driver

Cover screw (CS38)
connection

Healing Abutment



By hex driver



Healing abutment connection

Healing abutment (HAB402020L)
connection in thin gingiva

Surgical Kit Maintenance

Sterilization and Instrument Care Procedures


- Please follow legal regulations, as well as hygienic guidelines to prevent contamination and infection through prevention.
 - Please remember that you are responsible for the maintenance and sterility of your medical/dental products/device. It is important to use and follow proper cleaning, disinfection and sterilization procedures.
 - It is also important to follow the manufactures recommendation on the usage of drills. Please keep a log as to how many times the drills are used.
 - Drills are used per implant placed not per patient. Bone density determines the life of the drills.
 - Replace o-rings on adaptes and hex drivers, if worn and dried out.
 - Drills should be considered for replacement around 20-40 uses based on bone density.
- 01** All instruments immediately after use must be pre-soaked for a few minutes in a germicidal bath to loosen and prevent debris from attaching to instruments. Do not soak over-night.
 - 02** Scrub with a soft brush to remove any debris.
 - 03** For internal irrigation drills use a reamer or small gauge needle to cleanout drill internally.
 - 04** If using an ultrasonic cleaner, wrap drills in a 2 x 2 to prevent rubbing against each other.
 - 05** Rinse thoroughly under warm water.
 - 06** Clean all instrument trays with a germicidal cleaner prior to replacing instruments in kit.
 - 07** Dry completely and place back into kit.
 - 08** Always check for damage or corrosion after rinsing and drying.
 - 09** Seal the tray in a sterilization pouch.
 - 10** Sterilize using a steam autoclave in 121°C/250F for 30 minutes or refer to manufactures recommendations.
 - 11** Store in a dry area at room temperature.

Surgical Kit Maintenance


Maintenance Period for Surgical Drills

All surgical drills shall be replaced after approximately 40 uses based on bone density


● **Guide Drill**
(1000rpm, 30-45N · cm with Irrigation)




● **Final Drill**
(1000rpm, 30-45N · cm with Irrigation)




● **First Drill**
(1000rpm, 30-45N · cm with Irrigation)




● **Countersink Drill**
Depending on bone density, the rpm could be adjusted



● **Pilot Drill**
(1000rpm, 30-45N · cm with Irrigation)




● **Lindemann Drill (Guide, First)**
(1000rpm, 30-45N · cm with Irrigation)




Maintenance Period for DASK Drills

All DASK drills shall be replaced with new one after approximately 20 uses based on bone density


● **DASK Drill #1**
(800-1,200rpm, 30-45N · cm with Internal Irrigation)




● **DASK Drill #1**
(800-1,200rpm, 30-45N · cm with Internal Irrigation)




● **DASK Drill #2**
(800-1,200rpm, 30-45N · cm with Internal Irrigation)




● **DASK Drill #2**
(800-1,200rpm, 30-45N · cm with Internal Irrigation)



● **DASK Drill #3**
(800-1,200rpm, 30-45N · cm with Internal Irrigation)



● **DASK Drill #3**
(800-1,200rpm, 30-45N · cm with Internal Irrigation)



PROSTHETIC MANUAL

Prosthetic Introduction

Understanding the Implant and Prosthesis	12
Types of Abutments	13
Dual Abutment	14
Combi Abutment	15
Dual Milling / Anged / Temporay / Direct-Casting Abutmentt	16
Screw Abutment	17
Points to consider in Abutment Selection	18
Minimum Height Requirement for SuperLine Prosthetic Abutment	19

Prosthetic Procedure


Dual Abutment (Abutment Level Impression-Multiple Units)	20
Dual Abutment (Fixture Level Impression-Transfer Type, Multiple Units)	23
Dual Abutment (Fixture Level Impression -Pick Up Type, Multiple Units)	26
Combi Abutment (Abutment Level Impression-Multiple Units)	29
Dual Milling Abutment (Fixture Level Impression-Transfer Type, Single Units)	31
Screw Abutment (Abutment Level Impression-Transfer Type, Multiple Units)	33
Direct Casting Abutment (Fixture Level Impression-Single Units)	36
Angled Abutment (Fixture Level Impression-Pick Up Type, Single Units)	37
Temporary Abutment	38
Cementation Repair Method (SCRIP)	39
Ball Abutment	41
Dentium Prosthetics Instruction For Use	43
Dentium Implants Instruction For Use	45

Understanding the Implant and Prosthesis



Biological Connection

- The tapered conical hex connection between implant and abutment interface ensures hermetic sealing.
- The biological connection distributes the load to the fixture evenly. Therefore it may minimize bone loss.
- All implant diameters share the same internal connection. One abutment screw fits all abutments and fixtures.



Types of Abutment (Abutments are available in various diameters & gingival heights)

- Dual Abutment
 - Combi Abutment
 - Dual Abutment
 - Dual Milling Abutment
 - Angled Abutment (15° / 25°)
 - Direct-Casting Abutment
 - Temporary Abutment (Plastic & Titanium)
 - Screw Abutment
 - Ball Abutment
- Abutment level
Fixture level
Screw retained (Abutment level)
For denture use

Types of Abutment

One-Piece	Two-Pieces	
	 <p>Hex Non-hex</p>	 <p>Hex Non-hex</p>
 <p>Screw Abutment Abutment Level</p>	<p>[15°] [25°]</p>  <p>Hex Non-hex Hex Non-hex</p>	 <p>Hex Non-hex</p>
 <p>Cylinder</p>	<p>Titanium</p>  <p>Hex Non-hex</p>	<p>Plastic</p>  <p>Non-hex</p>
<p>Temporary Abutment Fixture Level</p>		

- Straight abutments are Dual and Combi.
- Depending on the insertion angle and position of the fixture, the Angled or Direct Casting abutment may be used.
- The Screw abutment can be used when prosthesis retrieval is anticipated.

Selection Guideline

Ideal emergence profile for each tooth



Dual Abutment



Dual Abutment

- It is possible to take an impression at both fixture level and abutment level.
(A dual abutment may be interchanged with a combi abutment.)
- For abutment level impressions, the same prosthetic procedures apply to both dual and combi abutments.
- For fixture level impressions, the abutment selection takes place on the master model.
- For fixture level impressions, a precise positioning jig for abutment may be required.
- Either hex or non-hex abutments may be used, according to operator's preference.

* If a cement retained restoration requires retrieval, cutting a hole in the occlusal surface would allow access to the screw to permit removal.

Hex / Non-hex


	Hex	Non-hex
Positioning Jig	Unnecessary	Required
Radiograph	Required	Unnecessary

Dual Abutment (Hex / Non-hex)

Diameter	G/H	Verticle angle
Ø4.5	1.0mm, 1.5mm, 2.5mm, 3.5mm, 4.5mm, 5.5mm	5°
Ø5.5	1.5mm, 2.5mm, 3.5mm, 4.5mm, 5.5mm	6°
Ø6.5	1.5mm, 2.5mm, 3.5mm, 4.5mm, 5.5mm	7°



Combi Abutment



- The Combi abutment is used when the implant position is optimal.
 - If the abutment selection is made in the mouth, gauge the thickness of mucosa with the depth gauge to measure the gingival height thus allowing the appropriate abutment height.
 - The Impression is taken with the snap cap
 - When using the Combi abutment, it remains in the mouth after impression. (DO NOT REMOVE OR CHANGE ITS POSITION).
 - Tighten abutment screw to 25 - 30 N·cm (retighten again before seating final prosthesis).
- * If the Combi abutment is too long it can be adjusted 1.5mm to the bottom of the laser mark on the vertical stack of the abutment. The Combi abutment has a short analog for the 1.5mm adjustment.
- * A resin jig can be made to record the reduction if reduced more than 1.5mm.

Combi Abutment Line Up

Diameter	G/H	Vertical angle
Ø4.5	1.0mm, 1.5mm, 2.5mm, 3.5mm, 4.5mm, 5.5mm	5°
Ø5.5	1.5mm, 2.5mm, 3.5mm, 4.5mm, 5.5mm	6°
Ø6.5	1.5mm, 2.5mm, 3.5mm, 4.5mm, 5.5mm	7°



Dual Milling / Angled / Temporary / Direct Casting Abutment




Dual Milling Abutment

- Impression is taken at fixture level.
- When using a non-hex abutment a precise seating jig should be used.
- Either hex or non-hex abutments may be used, according to operators preference.
- * if a cement retained restoration requires retrieval, cutting a hole in the occlusal surface would allow access to the screw for removal.

Angled Abutment

- The Angled Abutment is recommended when the restoration path of insertion is unfavorable in either anterior or posterior sites.
- Retention force can be increased through milling process.



Temporary Abutment

- Temporary abutments are available in titanium or plastic.
- The titanium abutment comes in hex and non-hex both with a gingival height of 1mm.
- The plastic abutment comes in diameters (4.5, 5.5, 6.5) with a gingival height of 3mm.

Direct-Casting Abutment

- Excellent for either single or for bridgework.
- Used as an esthetic custom made abutment.
- Used when angulation is not ideal and a standard abutment cannot be used.
- Used when there is inadequate inter-arch distance and a standard abutment cannot be used.
- A fixture level impression is taken, and the soft tissue contours can be supported.

Fixture Level Abutment (Hex / Non-hex)

Abutment	Diameter	G/H	Angle	
Dual Milling 	Ø4.0	1.5mm	×	
	Ø4.5	2.0mm		
	Ø5.5	2.5mm		
	Ø6.5	3.0mm		
Angled 	Ø4.5	2.0mm 4.0mm	15°/ 25°	
	Ø5.5	2.0mm 4.0mm	15°/ 25°	
	Plastic Temporary 	Ø4.5	3.0mm	×
		Ø5.5	3.0mm	
Ø6.5		3.0mm		
Direct-Casting Abutment 	Ø4.5	1.0mm	×	

Screw Abutment




If prosthesis repair is anticipated, use of a screw abutment retained prosthesis enables easy retrieval.

- Useful for connecting multiple units or if there is a preference for a screw retained prosthesis.
- Useful when respective long axes of implants differ. Each side tapers by 30° and this permits up to 60° divergence between two abutments.
- Useful if the prognosis of an adjacent restoration is not ideal thus permitting easy retrieval and modification of the restoration.

Ti-Retaining Screw (1.8mm - body diameter)

- Can minimize screw loosening due to increased approximal space.
 - Can endure various kinds of masticatory force.
- * 10N·cm of torque is recommended for the Ti-Retaining Screw



Screw Abutment

Diameter	G/H
Ø4.5	1.0mm, 1.5mm, 2.5mm, 3.5mm, 4.5mm, 5.5mm
Ø5.5	1.5mm, 2.5mm, 3.5mm, 4.5mm, 5.5mm



Points to Consider in Abutment Selection

Considerations in selecting an abutment

- Esthetic requirement
- Implant angulation
- Implant location
- Fixture installation depth (Gingival height)
- Interarch distance
- Prosthesis type
- Dentist & Dental technician's preference

Impression of Implant

According to the case the impression can be taken at abutment or fixture level.
Remove description and insert picture of transfer copings.

Fixture Level

1. Dual abutment
2. Dual milling abutment
3. Angled abutment (15°)
4. Direct casting abutment
5. Temporary abutment (Plastic & Titanium)


Abutment Level

1. Dual abutment
2. Combi abutment
3. Screw abutment

Abutment impression recommendation

Dual	Cementation type, screw-cementation type	Fixture level impression or abutment level impression
Combi	Cementation type	Abutment level impression
Angled	Cementation type, screw-cementation type	Fixture level impression
Screw	Screw retained type	Abutment level impression
Direct casting	Cementation type, screw-cementation type	Fixture level impression
Dual milling	Cementation type, screw-cementation type	Fixture level impression

Minimum Height Requirement for Prosthetic Abutment



※ Diagram above indicates the minimum height required for SuperLine/ IMPLANTIUM prosthetic abutment

Maximum amount of reduction allotted for adjustment

Combi Abutment

- Eliminate 3.0mm from the top level Combi abutment (laser marking:1.5mm)
Caution _ Damage may be caused to the screw if the abutment is reduced to less than 2.5mm above the gingival height.

Dual Abutment

- Preparation of the abutment top is possible as follows.

Gingival Height	Preparable Amount
1.5mm	2.0
2.5mm	3.0
3.5mm	4.0
4.5mm	5.0
5.5mm	6.0

Angled Abutment & Dual Milling Abutment

- Required minimum abutment height: at least 5.0mm above the Fixture top.

Direct-Casting Abutment

- Required minimum abutment height: at least 5.5mm above the Fixture top.

Screw Abutment

- The Screw abutment cannot be modified, however the casting abutment can be modified for interarch distance, taking reduction into consideration of the height of the retaining screw.

Combi Abutment

[Abutment Level Impression-Multiple Units]

Chairside



Remove Cover screw
[in case of second stage surgery]



Let soft tissue form around Healing
abutment



Select suitable Combi abutment,
then tighten it to 25-30N·cm.
Re-tighten after 15 minutes.



Select Impression
coping on Combi abutment.



Insert Impression coping over
Combi abutment firmly
[Snap-on Mechanism]



Inject impression material



Take impression



Impression coping comes off with
impression material



Fabricate provisional restoration, or
use Comport cap

Lab Side



Insert Lab analog into impression coping



Make sure Lab analog seats securely



Fabricate soft tissue model



Fabricate master cast



Seat burn-out cylinder securely into
Lab analog



Consider distance of opposing teeth,
Modify burn-out cylinder to its
proper height if needed

Combi Abutment

[Abutment Level Impression-Multiple Units]



Fabricate burn-out cylinder and plastic bar in preparation for wax-up



Wax up



Fabricate metal framework



Trim extended margin with rubber wheel



Metal framework and reamer



Eliminate the Lip remnant caused by 'snap-on' mechanism by reamer



Metal Framework after removal of the Lip remnant



Metal framework



Final prosthesis

Chairside



Seat final prosthesis and adjust occlusion

* If the Lab analog is trimmed due to limited inter-occlusal space in the lab, a reduction jig is necessary. Then slight modification of the abutment in the oral cavity may be necessary to reduce the height of the abutment

Dual Abutment

[Abutment Level Impression-Multiple Units]


Clinical Procedure



Chairside



Let soft tissue form around Healing abutment



Dual abutment (Hex/ Non-hex)




Select Dual abutment by diameter and gingival height




Tighten it to 25-30N ·cm. Re-tighten after 15 minutes



Insert Impression coping over abutment firmly [Snap-on Mechanism]



Inject impression material



Take impression



Impression coping comes off in the impression.



Fabricate provisional restoration, or use Comfort cap

Dual Abutment

[Abutment Level Impression-Multiple Units]

Laboratory Procedure



Lab Side



Insert Lab analog into impression coping



Make sure Lab analog seats securely into the impression coping [match flat side of both analog and coping]



Fabricate Soft tissue model



Fabricate master cast



Seat burn-out cylinder securely into Lab analog



Consider distance of opposing teeth, Modify burn-out cylinder to its proper height if needed



Fabricate burn-out cylinder and plastic bar in preparation for wax-up



Wax up



Fabricate metal framework

Dual Abutment

[Abutment Level Impression-Multiple Units]



Trim extended margin with rubber wheel



Metal framework and reamer



Eliminate the Lip remnant caused by 'snap-on' mechanism by reamer



Metal Framework after removal of the Lip remnant



Metal framework



Final prosthesis

SCRIP : Once an access hole has been created, it could be converted to a SCRIP (Screw & Cemented Retained Prosthesis).



Access hole is made when burn-out cylinder is used to do the wax up.



Appear extended margin around metal framework due to snap on mechanism.



Trim extended margin by rubber wheel



Metal framework and reamer



Eliminate the Lip remnant caused by 'snap-on' mechanism by reamer



Metal Framework after removal of the Lip remnant



Metal framework



Final prosthesis


Dual Abutment

[Fixture Level Impression-Transfer Type, Multiple Units]

Clinical Procedure



Chairside



Dual Abutment

[Fixture Level Impression-Transfer Type, Multiple Units]

Laboratory Procedure



Lab Side



Remove the impression coping from oral cavity and connect it with analog firmly



Attach the impression coping to the Alalog and insert into the impression firmly



Fabricate Soft tissue model



Fabricate master cast



Soft tissue condition after retrieval of impression coping



Measure gingival height with depth gauge



Select Dual abutment with proper diameter and gingival height



Verify the selected abutment by surveying (preparation is possible if necessary)



Fabricate positioning jig

Dual Abutment

[Fixture Level Impression-Transfer Type, Multiple Units]



Fabricate cap with pattern resin



Wax up



Metal framework



Final prosthesis

Chairside



Use the positioning jig to transfer the abutment in model to oral cavity. Then tighten it to 25-30N · cm. Re-tighten after 15 minutes.



Seat the final prosthesis and adjust occlusion. Place lab wax into opening site of abutment to protect screw head and then cement

SCRP-Lab Side



Make access hole in the resin cap by using a long impression coping transfer screw



Wax up



Metal framework



Final prosthesis

SCRP-Chairside



Use the positioning jig to transfer abutment in model to oral cavity then tighten it to 25-30N · cm. Re-tighten after 15 minutes.



Seat the final prosthesis and adjust occlusion. Place wax into opening site of the abutment prior to sealing with composite.

* In the process of seating the prosthesis, the prosthesis can be rebounded by gingival tissue. In that case it is advised to apply occlusal load on the prosthesis for 10-15 minutes.


Dual Abutment

[Fixture Level Impression -Pick Up Type, Multiple Units]

Clinical Procedure



Chairside



Dual Abutment

[Fixture Level Impression -Pick Up Type, Multiple Units]

Laboratory Procedure



Lab Side



Connect impression coping with analog firmly



Fabricate soft tissue model



Fabricate master cast



Select abutment with proper diameter and gingival height



Verify the selected abutment by surveying (preparation is possible if necessary)



Fabricate positioning jig



Fabricate cap with pattern resin



Wax up



Metal framework

Dual Abutment

[Fixture Level Impression -Pick Up Type, Multiple Units]

Chairside



Final prosthesis



Use positioning jig to transfer the abutment in model to oral cavity then tighten it to 25-30N · cm. Re-tighten after 15 minutes.



Seat the final prosthesis and adjust occlusion

*In the process of seating the prosthesis, the prosthesis can be rebounded by gingival tissue. In that case it is advised to apply occlusal load on the prosthesis for 10-15 minutes.

SCRP-Lab Side



Make access hole in the resin cap by using a long impression coping transfer screw



Wax up



Metal framework

SCRP-Chairside



Final prosthesis



Use positioning jig to transfer the abutment in model to the oral cavity then tighten it to 25-30N · cm. Re-tighten after 15 minutes.



Seat the final prosthesis and adjust occlusion

*In the process of seating the prosthesis, the prosthesis can be rebounded by gingival tissue. In that case it is advised to apply occlusal load on the prosthesis for 10-15 minutes.


Dual Milling Abutment

[Fixture Level Impression-Transfer Type, Single Unit]


Clinical Procedure



Chairside



Laboratory Procedure



Dual Milling Abutment

[Fixture Level Impression-Transfer Type, Single Unit]

Lab Side



Remove the impression coping from oral cavity and connect it with Lab analog firmly



Fabricate soft tissue model



Fabricate master cast



Select Dual milling abutment which has proper diameter



Abutment after milling process.



Fabricate positioning jig



Fabricate cap with pattern resin




Wax-up



Metal framework

Chairside



Final prosthesis



Use positioning jig to transfer the abutment in model to oral cavity then tighten it 25-30N • cm. Re-tighten after 15 minutes.




Seat final prosthesis and adjust occlusion

* In the process of seating the prosthesis, the prosthesis can be rebounded by gingival tissue. In that case it is advised to apply acclusal load on the prosthesis for 10-15 minutes.


Angled Abutment

[Fixture Level Impression-Pick Up Type, Single units]


Clinical Procedure



Chairside



Laboratory Procedure



Angled Abutment

[Fixture Level Impression-Pick Up Type, Single units]

Lab Side



Connect impression coping with analog firmly



Fabricate soft tissue model



Unscrew the impression coping screw, separate impression coping from the model



Fabricate master cast




Select and seat proper Angled abutment in master cast



Modify Angled abutment properly and fabricate positioning jig



Fabricate cap with pattern resin



Wax up



Metal framework

Chairside



Final prosthesis



Seat the Angled abutment using positioning jig



Insert final prosthesis and adjust occlusion


Direct Casting Abutment

[Fixture Level Impression-Single units]


Laboratory Procedure



Lab Side




Chairside




Temporary Abutment

[Fixture Level Impression-Pick Up Type, Single units]




Ti-Temporary Abutment Plastic Temporary Abutment

<Using Ti Cylinder>



Consider the opposing teeth height before seating the temporary abutment, Modify the abutment if needed and complete the temporary abutment prosthesis with direct resin


<Using Plastic Cylinder>




Screw Abutment

[Abutment Level Impression-Transfer Type, Multiple units]


Clinical Procedure




Cover Screw




Healing Abutment



Screw Abutment




Impression Coping
Transfer Type




Abutment Level
Impression


Chairside



Screw abutment with delivery holder



Select and seat appropriate screw abutment with delivery holder.



Tighten it to 25-30N·cm with S/B driver for screw abutment. Re-tighten after 15 minutes

Screw Abutment




Impression coping [transfer type] for Screw abutment



Seat impression coping on Screw abutment



Inject impression material



Take Impression




Inner-surface of impression



Seat comfort cap on the Screw abutment

Laboratory Procedure



Lab Analog Connection

Cylinder Connection

Modification

Crown Wax-Up

Final Restoration
Screw Retained

Lab Side



Remove the impression coping from oral cavity and connect it with analog firmly



Attach the impression coping to the Alalog and insert into the impression



Fabricate soft tissue model

Screw Abutment



Fabricate master cast



Remove impression coping



Connect the screw abutment cylinder then tighten it to 10N·cm with Ti-Retaining screw



Consider distance of opposing teeth, Modify cylinder to its proper height if needed



Fabricate burn-out cylinder and plastic bar in preparation for wax-up



Wax up



Fabricate metal framework



Eliminate the Lip remnant by reamer caused by 'snap-on' mechanism



Metal Framework after removal of the Lip remnant



Complete porcelain build up



Seat final prosthesis and adjust occlusion. Tighten it to 10 N·cm with Ti-Retaining screw

Cementation Repair Method (SCRMP)

In light of Implant Prosthesis:

- A screw type restoration helps simplify the prosthesis repair or insertion and removal of the prosthesis to any given situation.
- A dual abutment can be cement retained or screw retained.
- Combi abutment is only cement retained and occlusal hole is unnecessary.

In case of Screw Loosing or Prosthesis Repair is needed



In case of screw loosening and/or Prosthesis repair is needed



In order to unscrew, make access hole on the occlusal surface with bur



Unscrew, then remove the prosthesis from the oral cavity



Both cement retained prosthesis and abutment are removed



Finish the repair then re-seat into the oral cavity with new abutment screw



Tighten the prosthesis to 25-30N · cm with screw driver

* In case of screw abutment, Ti-Retain screw should be tighten to 10N · cm.



Fill the access hole with cotton



Then, fill the access hole with resin



Final prosthesis

Cementation Repair Method (SCRCP)

[Screw & Cement Retained Prosthesis]

Prosthesis separation from Abutment due to Cement Loss



Restore the separated prosthesis to the abutment in the oral cavity.



Unscrew the abutment screw to 25-30N·cm and remove prosthesis from the oral cavity completely.
* In case of screw abutment separation, Ti-Retain screw should be unscrewed to 10N·cm.



Apply cement to the prosthesis



Adhere the prosthesis to abutment and clean out remnant cement, Fill the access hole with cotton and resin

In case of prosthesis loosening and adding to the interproximal surface is necessary



Adding to the inter proximal contact due to loosening



Make access hole using bur



Unscrew then remove the cemented prosthesis with abutment in the oral cavity



Add resin to the contact if needed



Insert the prosthesis in the oral cavity and screw it in. Perform light curing, then polish the contact area.



Replace the prosthesis into the oral cavity and tighten with a new abutment screw, fill access hole with cotton and resin.



Final Prosthesis after repair

Ball Abutment



Connect Ball abutment with fixture



Seat impression coping into Ball abutment



Make individual tray for denture impression



Inject impression material



Take impression with individual tray



Remove the tray from the oral cavity



Inner surface of the impression material




Ball analog



Insert analog into the impression material securely



After analog insertion



Fabricate master cast




Fabricate denture with common method



Connect female socket with analog



Reduce denture inner surface to place the female socket



Examine the interference between inner surface of ball analog and female socket

Ball Abutment




Apply the resin with brush into the hole



Remove the denture from the model after initial setting of resin



Add resin with brush around the female socket



After the resin sets, trim the remnant resin from the denture.

Dentium Prosthetics Instruction For Use

Product Description

Dentium Prosthetics are devices made of titanium alloy intended for use as an aid in prosthetic rehabilitation. They consist of Angled Abutments, Ball Abutments, Combi Abutments, Direct-Casting Abutments, Dual Abutments, Dual Milling Abutments, Healing Abutments, Screw Abutments, Temporary Abutments and Abutment Screws. Most abutments are partially TiN coated. Abutments are supplied non-sterile and are intended to be sterilized prior to use (see **Sterilization**).

Indication for Use

Dentium Prosthetics are intended for use as an aid in prosthetic rehabilitation.

Contraindications

- Blood dyscrasia not under control
- Alcoholic, drug use
- Poor osseous patients unable to implant
- Other periodontal and surgical operation prohibited patients
- Metal allergy

Warnings

- Implant surgery and restoration involves complex dental procedures which require specialized training. Training is strongly recommended prior to prosthetic use.
- Improper technique and/or improper patient selection can result in implant failure and/or loss of bone around the implant site.

Precautions

Thorough screening of implant candidates is critical to the success of the implant process. Appropriate radiographic examination should be utilized to determine adequacy of bone, periodontal status, and the location of important anatomical landmarks. Post-implantation radiographs are also required to determine the progression of osseointegration.

Exposure to magnetic resonance imaging, radiation, and chemotherapy may impact the health of the implant. Dental implant patients should be instructed to consult with their physicians prior to undergoing such treatment options. The use of electrosurgical instruments or lasers around the titanium implants and their abutments is not recommended due to the risk of electrical shock and/or burns.

Sterilization

Dentium Prosthetics are supplied non-sterile. Open the pouch and take out the vial. Pull out the cap and remove product from the vial. Prior to use, sterilize by steam heat 121°C (250°F) for 30 minutes.

Instructions for Use (Prosthetic Procedure)

Dentium Prosthetics are part of Dentium's Dental Implant System and are only intended for use with Dentium Implants and Surgical Instruments.

Dentium Prosthetics Instruction For Use

Caution: Refer to the Surgical and Prosthetic Manual for detailed instructions regarding prosthetic procedures.

Warning: Evidence of osseointegration should be confirmed radiographically prior to loading a fixture to avoid movement of the fixture and possible loss of the implant.

Warning: Do not overload implants. Excessive force may cause fracture of the implant/abutment. Careful treatment planning is required to ensure a sufficient number of implants are used to create adequate contact surface and avoid overloading.

Adverse Effects

The implant surgery has known risks, including localized swelling, separation of sutured area, ulceration, bleeding, infection, and localized tenderness of short duration. In some cases, there can be numbness in the area of the surgery which is typically temporary but may be permanent if there is nerve damage. Most of these risks are managed by standard pre and post operative care prescribed by the dentist or surgeon. Patients should be instructed to seek the services of a trained dental professional if there are any signs of a problem with the implant (such as looseness, infection, or exudates around the implant, pain or other symptoms not expected).

Loss of the implant is possible due to a lack of osseointegration for the following reasons: quality or quantity of bone is insufficient, placement or alignment of the implant is not optimal, patient oral hygiene, health or habits are not optimal, or the implant is loaded prior to adequate osseointegration. As with all implant systems, there is the risk of bone loss around the implant over time which can result in revision or removal of the implant.

Cautions

- Federal law restricts this device to sale by or on the order of a dentist or physician.
- Do not use the product if the packaging is damaged or opened.
- Do not use beyond the stated expiration date.
- Product is designed for single use only and is not to be reused. Opened, unused product may not be returned to the company.

Delivery, Storage and Handling

- Prosthetic: NON-STERILE
 - Store product at room temperature.
-

Dentium Implants Instruction For Use

Product Description

IMPLANTIUM® and SuperLine™ fixtures are titanium dental implants with Sandblasting with Large grit and Acid etching surface treatment designed to enhance osseointegration, a firm connection between the bone and the surface of the implant, to provide a stable support for prosthetic devices.

Indication for Use

IMPLANTIUM® and SuperLine™ fixtures are intended for use in dental implant surgery. The fixtures are surgically placed in the upper or lower jawbone to provide support for dental prostheses including single tooth replacements, bridgework or complete dentures.

Contraindications

IMPLANTIUM® and SuperLine™ fixtures should not be used in patients with insufficient bone width or height. Implant failure can occur in patients with insufficient bone, poor bone quality, poor oral hygiene, habits such as heavy smoking, tobacco abuse, alcohol or drug abuse, or generalized disease such as blood disorders or diabetes. Fixtures should not be used in patients who are allergic to titanium or titanium alloys.

Warnings

Implant surgery and restoration involve complex dental procedures which require specialized training. Training is strongly recommended prior to implant use. Improper technique and/or improper patient selection can result in implant failure and/or loss of bone around the implant site.

Precautions

Thorough screening of implant candidates is critical to the success of the implant process. Appropriate radiographic examination should be utilized to determine adequacy of bone, periodontal status, and the location of important anatomical landmarks. Post-implantation radiographs are also required to determine the progression of osseointegration.

Exposure to magnetic resonance imaging, radiation, and chemotherapy may impact the health of the implant. Dental implant patients should be instructed to consult with their physicians prior to undergoing such treatment options. The use of electrosurgical instruments or lasers around the titanium implants and their abutments is not recommended due to the risk of electrical shock and/or burns.

Adverse Effects

The implant surgery has known risks, including localized swelling, separation of sutured area, ulceration, bleeding, infection, and localized tenderness of short duration. In some cases, there can be numbness in the area of the surgery which is typically temporary but may be permanent if there is nerve damage. Most of these risks are managed by standard pre and post operative care prescribed by the dentist or surgeon. Patients should be instructed to seek the services of a trained dental professional if there are any signs of a problem with the implant (such as looseness, infection, or exudates around the implant, pain or other symptoms not expected).

Loss of the implant is possible due to a lack of osseointegration for the following reasons: quality or quantity of bone is insufficient, placement or alignment of the implant is not optimal, patient oral hygiene, health or habits are not optimal, or the implant is loaded prior to adequate osseointegration. As with all implant systems, there is the risk of bone loss around the implant over time which can result in revision or removal of the implant.

Dentium Implant Instruction For Use

Cautions

- Federal law restricts this device to sale by or on the order of a dentist or physician.
- Open fixture vial with cap up to ensure the fixture remains in place to attach to the Implant Driver.
- Refer to the Surgical and Prosthetic Manual for detailed instructions regarding the surgical procedure.
- Product is gamma sterilized and will remain sterile as long as the packaging is unopened and undamaged. Do not use the product if the packaging is damaged or opened. Do not use beyond the stated expiration date. Product is designed for single use only and is not to be resterilized. Opened, unused product may not be returned to the company.

Delivery, Storage and Handling

- Outer Box: NON-STERILE
- Blister Packaging: NON-STERILE
- Outer Vial: STERILE-R
- Inner Vial: STERILE-R
- Implant: STERILE-R

Store product at room temperature. Open product using aseptic technique to avoid contamination after correct size has been determined prior to implantation. Handle the fixture carefully to avoid damaging the fixture surface.

Additional product labels containing the product LOT number and REF number are provided with the fixture. These labels are to be placed in the patient's chart and medical records to ensure complete traceability of the implant for future reference or questions.

To obtain additional information on IMPLANTIUM® and SuperLine™ products, please contact:

Dentium USA
Cypress, CA 90630
Tel: 1-877-304-6752
www.dentiumusa.com

DentiumUSA
Developed by Clinicians for Clinicians

SuperLine[™] & IMPLANT[®]

DentiumUSA

Developed by Clinicians for Clinicians

Specifications are subject to change without prior notice.
Some products to be launched in the market after necessary approvals are also listed in this catalog.

Dentium USA

6761 Katella Ave. Cypress, CA 90630 Tel. 1-877-304-6752 Fax. 714-226-0019

Website / E-mail

www.dentiumusa.com / info@dentiumusa.com

PN2433-02 03/10