

Dentium



More Speed

IMPLANTIUM

Clinical Success & Simplicity

1. Optimal fixation screw

- Minimal bone loss due to adequate stimulation to marginal cortical bone.
- Good marginal bone remodelling with minimal bone loss.

2. Biological connection

- Conical connection between fixture and abutment interface ensures hermetic sealing.
- Distributes the load to the surrounding cortical bone evenly.

3. Surface

- S.L.A. surface promotes early osseointegration by fast osteoblast attachment.
- Machined surface top bevel for healthy gingiva & simple 2nd stage surgery.

4. Abutment partial Ti-N coating

- Ti-N coating leads to better gingiva aesthetics through gingival translucency at abutment interface.



Australia T +61-2-9874-0183	Lebanon T +961-1254515	Saudi Arabia T +961-3-10-1211
China T +86-10-6280-0505	Mexico T +52-664-638-3398	Singapore T +65-6767-1311
Hong Kong T +852-3741-1222	Poland T +48-22-313-0885	Taiwan T +886-2-2706-1279
Iran T +98-21-4426-5178	Romania T +40-21-410-8326	Uzbekistan T +99871-144-3055
Korea T +82-2-501-8560	Russia T +7(495)-232-3952	

IMPLANTIUM

www.implantium.com

Dentium USA 11075 Knott Ave., Suite A, Cypress, CA 90630, USA T +1-714-895-7278 F +1-714-890-6360